

## Vegetables



# Reagents and Standards for Vegetables



As per Official Methods of Analysis of AOAC International, 19th Edition, 2012 Volume 2.  
Dr George W. Latimer, Jr.

## AOAC Official Method 974.24 Oxalic Acid in Canned Vegetables

Chapter	Product No.	Description	Pack Size
42.1.18A.b	PP2001F	Analytical Volumetric Solution Potassium Permanganate 0.05N, 0.01M	1L
42.1.18A.c	DBC06-960	Dissolution Media Acetate Buffer pH 4.5 Concentrate (Dilute to 25L) (6x961.5ml)	6X 1L

## AOAC Official Method 956.06 Aldehydes as Acetaldehyde in Frozen Vegetables

Chapter	Product No.	Description	Pack Size
42.3.06B.a	T20051	Analytical Volumetric Solution Sodium Thiosulphate 0.05N, 0.05M	1L
42.3.06B.b	I2005F	Analytical Volumetric Solution Iodine 0.1N, 0.05M	1L

## AOAC Official Method 947.11 Catalase in Frozen Vegetables

Chapter	Product No.	Description	Pack Size
42.3.07A.d	3002900	European Pharmacopoeia Reagent 0.01M Iodine	1L

## AOAC Official Method 963.27 Peroxidase in Frozen Vegetables

Chapter	Product No.	Description	Pack Size
42.3.08A.a	S21001	Analytical Volumetric Solution Sodium Hydroxide 1.0N, 1.0M	1L
42.3.08A.c	NACL2WV5	Sodium Chloride 2%	500ml
42.3.08A.d	SU2200J	Analytical Volumetric Solution Sulphuric Acid 4.0N, 2.0M	2.5L
42.3.08A.d	KI2101	Analytical Volumetric Solution Potassium Iodide 1.0N, 1.0M	1L
42.3.08A.d	T20101	Analytical Volumetric Solution Sodium Thiosulphate 0.1N, 0.1M	1L
42.3.08A.e	T20011	Analytical Volumetric Solution Sodium Thiosulphate 0.01N, 0.01M	1L

# Reagecon's Conductivity Standards

Product No.	Description	Pack Size
CSKC13	1.30µS/cm @ 25°C	250ml
CSKC136	1.30µS/cm @ 25°C	6x250ml
CSKC3	3µS/cm @ 25°C	250ml
CSKC5	5µS/cm @ 25°C	500ml
CSKC10	10µS/cm @ 25°C	500ml
CSKC20	20µS/cm @ 25°C	500ml
CSKC238	23.8µS/cm @ 25°C	500ml
CSKC25	25µS/cm @ 25°C	500ml
CSKC50	50µS/cm @ 25°C	500ml
CSKC84	84µS/cm @ 25°C	500ml
CSKC100	100µS/cm @ 25°C	500ml
CSKCS	147µS/cm @ 25°C	500ml
CSKC150	150µS/cm @ 25°C	500ml
CSKC185	185µS/cm @ 25°C	500ml
CSKC200	200µS/cm @ 25°C	500ml
CSKC250	250µS/cm @ 25°C	500ml
CSKC300	300µS/cm @ 25°C	500ml
CSKC400	400µS/cm @ 25°C	500ml
CSKC500	500µS/cm @ 25°C	500ml
CSKC718	718µS/cm @ 25°C	500ml
CSKC1000	1000µS/cm @ 25°C	500ml
CSKCL	1413µS/cm @ 25°C	500ml
CSKC2M	2000µS/cm @ 25°C	500ml
CSKC2500	2500µS/cm @ 25°C	500ml
CSKC3M	3000µS/cm @ 25°C	500ml
CSKC5M	5,000µS/cm @ 25°C	500ml
CSKC7M	7,000µS/cm @ 25°C	500ml
CSKC10M	10,000µS/cm @ 25°C	500ml
CSKC12880	12,880µS/cm @ 25°C	500ml
CSKC20M	20,000µS/cm @ 25°C	500ml
CSKC30M	30,000µS/cm @ 25°C	500ml
CSKC40M	40,000µS/cm @ 25°C	500ml
CSKC50M	50,000µS/cm @ 25°C	500ml
CSKC60M	60,000µS/cm @ 25°C	500ml
CSKC80M	80,000µS/cm @ 25°C	500ml
CSKC100M	100,000µS/cm @ 25°C	500ml
CSKC111800	111,800µS/cm @ 25°C	500ml
CSKC150M	150,000µS/cm @ 25°C	500ml
CSKC200M	200,000µS/cm @ 25°C	500ml
CSKC300M	300,000µS/cm @ 25°C	500ml
CSKC350M	350,000µS/cm @ 25°C	500ml
CSKC400M	400,000µS/cm @ 25°C	500ml
CSKC450M	450,000µS/cm @ 25°C	500ml
CSKC500M	500,000µS/cm @ 25°C	500ml

## Pesticide Standard's in Vegetables

Product No.	Description	Pack Size
REPET180	Pesticide Single Component Standard 1000ug/ml Chlorpyrifos	1ml
REPET148	Pesticide Single Component Standard 1000ug/ml Lindane (HCHgamma)	1ml